







VistaLab Catalog #	Color	Volume Range	Dispensing Increments	Accuracy	Precision
1065-0020 1065-0020L*		0.5 – 20 μ L	0.05 μ L	\pm 1.0% at 20 μ L \pm 1.5% at 10 μ L \pm 3.0% at 2 μ L	0.30% at 20 μ L 0.50% at 10 μ L 1.80% at 2 μ L
1065-0125 1065-0125L*		2 – 125 μ L	0.5 μ L	\pm 0.8% at 125 μ L \pm 0.8% at 62.5 μ L \pm 2.9% at 12.5 μ L	0.15% at 125 μ L 0.24% at 62.5 μ L 1.00% at 12.5 μ L
1065-0250 1065-0250L*		5 – 250 μ L	1 μ L	\pm 0.8% at 250 μ L \pm 0.8% at 125 μ L \pm 2.5% at 25 μ L	0.15% at 250 μ L 0.25% at 125 μ L 1.00% at 25 μ L
1065-1250 1065-1250L*		25 – 1250 μ L	5 μ L	\pm 0.75% at 1250 μ L \pm 0.8% at 625 μ L \pm 2.8% at 125 μ L	0.15% at 1250 μ L 0.20% at 625 μ L 0.60% at 125 μ L
1065-05mL 1065-05mLL*		100 – 5000 μ L	10 μ L	\pm 0.6% at 5000 μ L \pm 0.6% at 2500 μ L \pm 2.4% at 500 μ L	0.16% at 5000 μ L 0.20% at 2500 μ L 0.60% at 500 μ L
1065-10mL 1065-10mLL*		0.2 – 10mL	10 μ L	\pm 0.6% at 10mL \pm 1.0% at 5mL \pm 2.4% at 1mL	0.16% at 10mL 0.20% at 5mL 0.30% at 1mL

