

Fixed Volume Models

VistaLab Catalog #		Color	Volume	Accuracy	Precision
Single Stroke	Two Stroke				
	1257-0001	●	1 μ L	\pm 2.7%	1.80%
	1257-0002	●	2 μ L	\pm 1.5%	0.90%
1157-0005	1257-0005	●	5 μ L	\pm 1.5%	0.60%
1157-0010	1257-0010	●	10 μ L	\pm 1.5%	0.50%
1157-0020	1257-0020	●	20 μ L	\pm 1.0%	0.30%
1157-0025	1257-0025	●	25 μ L	\pm 1.4%	0.60%
1157-0050	1257-0050	●	50 μ L	\pm 0.8%	0.24%
1157-0100	1257-0100	●	100 μ L	\pm 0.8%	0.15%
1157-0200	1257-0200	●	200 μ L	\pm 0.8%	0.15%
1157-0250	1257-0250	●	250 μ L	\pm 1.2%	0.30%
1157-0300	1257-0300	●	300 μ L	\pm 1.0%	0.25%
1157-0500	1257-0500	●	500 μ L	\pm 0.8%	0.20%
1157-1000	1257-1000	●	1000 μ L	\pm 0.8%	0.15%
1157-0999	Ovation Fixed Volume Starter Pack (Includes one 10 μ L, one 100 μ L, one 200 μ L, one 1000 μ L, an Ovation Pipette Stand, and 3 racks of tips.)				



Custom volume settings are available upon request