

Tips with a Perfect Fit.

4060-1002 Catalog number

Use this tip chart to determine which tips are most appropriate for use with your Ovation and MLA pipettes.

		0.5-20µL	2-125µL 5-250µL	25-850µL 25-1250µL		
Ovation – electronic single channel and multichannel models:						
Ovation – electronic macro models:					100-5000µL	0.2-10mL
Ovation – adjustable volume models:		0.2-2µL 1.0-10µL 2.0-20µL	10-100µL 20-200µL	100-1000µL		
Ovation – fixed volume models:		1µL to 20µL	21µL to 200µL	201µL to 1000µL		
MLA models:			2µL to 200µL	201µL to 1000µL	1mL to 5mL	6mL to 10mL
TIP SIZE:		MICRO	SMALL	LARGE	5mL	10mL
NON-STERILE	VistaRak 192 tips/rack, 5 racks	4060-1002	4060-2004	4060-3004		
	VistaStak, 192 tips/layer, 5 layers (small size) or 3 layers (micro and large size)	4060-9024	4060-9025	4060-9026		
	Stacked Rack, 200 tips/layer, 5 layers (small size) or 3 layers (large size)		9025 ¹	9026 ¹		
	Stacked Rack, Trace Metal Certified, 200 tips/layer 5 layers (small size) or 3 layers (large size)		9022 ¹	9023 ¹		
	Box 100 tips/box (5mL) or 35 tips/box (10mL)				9048	9050
	VistaClear Box 35 tips/box					4058-6100
	VistaBulk™, 1000 tips/bag, or 250 tips/bag (5mL) or 100 tips/bag (10mL)	4058-1000	4058-2000	4058-3000	4058-5000	4058-6000
	Protectainer™ Bulk Pack, 1000 tips (small) or 750 tips (large)		4025 ²	4026 ³		
	Econo-Pak™ Bulk Pack, 1000 tips		4225 ²	4226 ³		
STERILE	VistaRak, Sterile, Pyrogen-free, RNase/DNase certified, 192 tips/rack, 5 racks	4060-1032	4060-2132	4060-3132		
	VistaClear Box – Sterile, 60 tips/box (5mL) or 35 tips/box (10mL)				4058-5102	4058-6102
	VistaTip™ Individually Wrapped Sterile, 200 tips		2025 ²	2026 ³		
	VistaTip Individually Wrapped, Sterile, Pyrogen-free, 200 tips		2027 ²	2028 ³		
FILTERED	VistaRak, Filtered, Sterile, Pyrogen-free, RNase/DNase certified, 192 tips/rack, 5 racks † 2µL tip for Ovation 0.2-20µL, Ovation 1µL and Ovation 2µL only	4060-1332 4060-1333†	4060-2332	4060-3332		
	VistaClear Box – Filtered, Sterile, Pyrogen-free, RNase/DNase certified, 60 tips/box (5mL) or 35 tips/box (10mL)				4058-5133	4058-6133
	VistaTip Individually Wrapped, Filtered, Sterile, Pyrogen-free, RNase/DNase certified, 50 tips				4058-5332	4058-6332



VistaRak™



VistaStak™



Stacked Rack



VistaClear™ Box

- Notes:**
1. Not for use with Ovation Multichannel Pipettes
 2. Not for use with 5–250µL Electronic Single Channel model
 3. Not for use with 25-1250µL Electronic Single Channel model
 4. Filtered tips contain a unique hydrophobic filter which acts as a barrier to aqueous liquids and aerosols, protecting the pipette and sample from trace amounts of carryover.
 5. All 5mL tips are graduated, marked at 0.1mL increments